



Blood Culture Identification (BCID) Verification Panel (Inactivated) - Catalog #8201*

Contains three vials of inactivated, pooled microorganism material and three vials of hydration buffer. Each vial contains enough material for five uses and can be used for up to 72 hours when refrigerated. For use with the The BioFire® FilmArray® Blood Culture Identification (BCID) Panel.

Blood Culture Identification (BCID) Control Panel (Inactivated) - Catalog #8215*

Contains six inactivated microorganism pellets, each pooled with 27 targets. For use with the The BioFire® FilmArray® Blood Culture Identification (BCID) Panel.

Respiratory (21 Targets) Control Panel - Catalog #8217*

Contains six inactivated microorganism pellets, each pooled with 21 targets. For use with The BioFire® FilmArray® Respiratory Panel 2 (RP2).



Accurate, Reliable Results

- Each panel contains the recommended strains for a particular assay in an easy-to-use format
- Designed with known performance on the BioFire® FilmArray® Systems
- Individually packaged and pre-pooled to simplify use and avoid cross contamination
- · Convenient test ready format saves you time and money
- · One catalog number for fast, easy ordering
- For in vitro diagnostic use (IVD)
- · Room temperature storage
- Technical Support experts available for guidance

*Live cultures also avaiable. Each set contains the recommended strains in the KWIK-STIK™ format. Learn more at microbiologics.com.

